

# Select Your Bullet Blender

## Start with sample size

*Choose the Bullet Blender that's right for your lab*



**up to 0.3g or 0.3 mL**

1.5 mL tubes  
24 samples



**Storm Pro**  
small footprint



**up to 1 g or 1 mL**

5 mL tubes  
12 samples



**5E Pro**  
small footprint



**up to 3 g or 3 mL**

50 mL tubes  
8 samples



**50 Gold**  
Maintains samples at 4°C  
Very quiet operation



**Gold**

Maintains samples at 4°C  
Very quiet operation



**5E Gold**

Maintains samples at 4°C  
Very quiet operation



Need a protocol? Email [tech@midsci.com](mailto:tech@midsci.com)