



TECHNICAL DATA SHEET

DEEP WELL PLATE



MID-KF-COMB-96

GENERAL SPECIFICATIONS

NUMBER OF WELLS:	96
MATERIAL:	Polypropylene
COLOR:	Natural
TEMPERATURE RANGE:	-80°C to 120°C
AUTOCLAVABLE:	121C at 15 psi, 15 min.
AUTOMATION COMPATIBLE:	Yes

CRITICAL DIMENSIONS

OVERALL HEIGHT:	45.00mm
TOP (L x W):	127.30mm x 85.30mm
WALL THICKNESS:	0.66mm
WELL A1 (X,Y):	14.15mm, 11.15mm
SINGLE WELL DEPTH:	44.11mm
WELL DIA. (TOP):	6.51mm
WELL DIA. (BOT):	4.25mm

PACKAGING CONFIGURATIONS

PACKAGING QUANTITY (CS):	100; 2 per Pack
CASE DIMENSIONS (IN.):	10.7 x 10.75 x 17.75
CASE WEIGHT (LBS.):	13.3

STORAGE

Keep out of sunlight. Store in a cool dry place. Keep away from sources of ozone or ignition.



COMPATIBLE PLATE OPTIONS

ITEM NUMBER:	DESCRIPTION:
MID-KF22US-9-N	2.2mL 96-Square well with Round V Bottom
MID-KF05US-9-N	200ul 96-Square well with Round V Bottom (elution)

All plates can be sealed with adhesive or heat seal films.