

**SPEC SHEET ABT-115V-1786**

<b>Product Name:</b>	ABT-115V-1786
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Ultra Low Temperature Freezer</li> <li>• 17 Cu.Ft. Capacity</li> <li>• One solid swing exterior door, right hinged, not reversible</li> <li>• Five insulated inner doors</li> <li>• Microprocessor controller specifications</li> <li>• Touch pad keys and four line display for easy programming</li> <li>• Large LED temperature readout</li> <li>• Audible and visual alarms with silence key for audible alarms</li> <li>• Programmable high and low temperature alarm set points</li> <li>• Programmable access code to limit access to set points and controller configuration</li> <li>• Automatic voltage boost to compensate for low line voltages</li> <li>• Normally open, normally closed alarm contacts and 4-20mA signal available to connect to central alarm system</li> <li>• RS-232 serial output</li> <li>• Battery backup on electronics during power outages</li> <li>• Temperature range: -50°C to -86°C</li> <li>• Four stainless steel shelves, adjustable on ½" centers</li> <li>• Cam action handle for one hand latching of exterior door.</li> <li>• Handle supplied with key lock and door switch for audible and visual alarms</li> <li>• Four heavy duty swivel casters with locks</li> <li>• Vacuum relief port for quick access after door closures</li> <li>• CFC-free refrigerant</li> <li>• Two year parts and labor warranty, plus an additional three year compressor parts warranty</li> <li>• Exterior dimensions: 37.6"W x 36.3"D x 79.6"H</li> <li>• Interior dimensions: 23.4"W x 23"D x 50"H</li> <li>• 2" box capacity: 320 ea.</li> <li>• Upright freezer rack capacity: 20 ea. UF442 racks</li> <li>• 115V, 60Hz, 16 amp, HP: ½ x ½</li> <li>• UL/C-UL Listed</li> <li>• NEMA plug configuration: NEMA 5-20</li> <li>• Cascade refrigeration system</li> <li>• Shipping weight: 700 lbs.</li> </ul>