

SPEC SHEET **ABT-HC-CP-16**

Product Name:	ABT-HC-CP-16
Features:	<ul style="list-style-type: none"> • Premier Glass Door Chromatography Refrigerator • 16 Cu. Ft. capacity • One swing glass door, self-closing, right hinged, not reversible • Microprocessor Temperature Controller • Adjustable operating temperature range: 1°-10°C • Digital Temperature Display and Alarm Module • Sample simulated digital temperature display; .1° Resolution • Adjustable temperature display in Celsius or Fahrenheit • Audible and visual high/low temperature alarms • Sensor error alarm • Power failure alarm • Resettable 1-point Min/Max History • Remote alarm contacts • 2 probes; 1 in air and 1 in sample bottle • Alarm validation (audible, visual, and remote alarm) • Alarm mute with ring back • Three adjustable shelves (adjustable in ½" increments) • Forced draft air circulation • Keyed door lock • One GFCI protected duplex electrical outlet • Two 2" access ports with spring-loaded, self-closing gasketed metal covers • Magnetic door gasket for positive seal • Heavy duty door handle • LED interior lights are safety shielded and switch controlled • Swivel Casters (3 ½") Probe access port (3/4") • Advanced hydrocarbon detector

SPEC SHEET **ABT-HC-CP-16**

Product Name:	ABT-HC-CP-16
Features:	<ul style="list-style-type: none"> • White powder-coated steel interior and exterior • High density urethane foam cabinet insulation • Hydrocarbon, natural refrigerant (R290) • Two year parts and labor warranty, plus an additional three compressor parts warranty • Exterior dimensions: 26 1/2" W* x 29 3/4" D x 79" H with casters; 75 1/2" H without casters. *Width w/access ports closed; 34 1/2" width w/access ports open • Cycle defrost • 115V, 60 Hz, 3 Amps, 1/5 HPETL listed • Shipping weight: 283 lbs. • Note: This unit requires two (2) outlets