

SPEC SHEET **ABT-HC-CP-36**

Product Name:	ABT-HC-CP-36
Features:	<ul style="list-style-type: none"> • Premier Glass Door Chromatography Refrigerator • 36 Cu. Ft. capacity • Two swing glass doors, self-closing, right hinged, not reversible • Microprocessor Temperature Controller • Adjustable operating temperature range: 1°-10°C • Digital Temperature Display and Alarm Module • Sample simulated digital temperature display; .1° Resolution • Adjustable temperature display in Celsius or Fahrenheit • Audible and visual high/low temperature alarms • 2 probes; 1 in air and 1 in sample bottle • Sensor error alarm, Power failure alarm • Resettable 1-point Min/Max History • Remote alarm contacts • Alarm validation (audible, visual, and remote alarm) • Alarm mute with ring back • Eight adjustable shelves (adjustable in 1/2" increments) • Keyed door lock • One GFCI protected duplex electrical outlet • Two 2" access ports with spring-loaded, self-closing gasketed metal covers • Forced draft air circulation • Magnetic door gasket for positive seal, Heavy duty door handle • LED interior lights are safety shielded and switch controlled • Swivel Casters (3 1/2"), Probe access port (3/4")

SPEC SHEET ABT-HC-CP-36

Product Name:	ABT-HC-CP-36
Features:	<ul style="list-style-type: none"> • Advanced hydrocarbon detector • White powder-coated steel exterior and interior • High density urethane foam cabinet insulation • Hydrocarbon, natural refrigerant (R290) • Two year parts and labor warranty, plus an additional three compressor parts warranty • Exterior dimensions: 41 1/8" W* x 34 3/4" D x 81 3/4" H with casters; 78 1/4" H without casters. *Width w/access ports closed; 49 1/8" width w/access ports open • Cycle defrost • 115V, 60 Hz, 4.5 Amps, 1/4 HP • ETL listed • Shipping weight: 405 lbs. Note: This unit requires two (2) outlets