

SPEC SHEET ABT-HC-EFP-14

Product Name:	ABT-HC-EFP-14
Features:	<ul style="list-style-type: none"> • Standard Hazardous Location (Explosion Proof) Freezer • 14 Cu. Ft. Capacity • One solid swing door, right hinged, not reversible • Microprocessor Temperature Controller • Adjustable operating temperature range: -15°C – -25°C • Six fixed shelves • Keyed door lock • Magnetic door gasket for positive seal • Heavy duty door handle • Leveling legs • White textured cabinet and door • High density urethane foam cabinet and door insulation • Hydrocarbon, natural refrigerant (R290) • Exterior dimensions: 23 5/8"W x 24 3/8"D x 73"H • Manual defrost • Not equipped with line cord and plug; must be hardwired in conduit to 120 VAC power source • 115V, 60Hz, 3.5 Amps, 1/3 HP • ETL listed • Two year parts and labor warranty • Meet NFPA and OSHA guidelines for hazardous locations, Class 1, Division II, Group C & D. These units are engineered for use in environments where volatile / explosive conditions could potentially exist. For storage of flammable or volatile materials within typical laboratory or non-volatile environments, see Flammable Material Storage units • Shipping weight: 195 lbs.