

**SPEC SHEET ABT-HC-SSP-49FA**

<b>Product Name:</b>	ABT-HC-SSP-49FA
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Premier Series Stainless Steel Freezer</li> <li>• 49 Cu. Ft. capacity</li> <li>• Two swing solid doors, self-closing, not reversible</li> <li>• 6 ea. Adjustable Shelves</li> <li>• Digital display microprocessor temperature controller</li> <li>• Adjustable operating temperature: -15°C to -25°C</li> <li>• Audible and visual Hi/Lo temperature alarm</li> <li>• 2 probes; 1 in air and 1 in sample bottle</li> <li>• Probe attached to evaporator for defrost termination sensing</li> <li>• Sensor error alarm</li> <li>• Continuous product (simulated) temperature display</li> <li>• Remote alarm contacts</li> <li>• Forced draft air circulation</li> <li>• Keyed door lock</li> <li>• Magnetic door gaskets for positive seal</li> <li>• Integrated door handle</li> <li>• LED interior lights are safety shielded and switch controlled</li> <li>• Probe access port (3/4")</li> <li>• Stainless steel exterior</li> <li>• White-powder coated aluminum interior with 304 stainless steel floor</li> <li>• Easily cleaned smooth interior</li> <li>• Hydrocarbon natural refrigerant (R290)</li> <li>• 2 year parts and labor warranty, plus an additional 3 year compressor parts warranty</li> <li>• Auto defrost</li> </ul>

**SPEC SHEET**    **ABT-HC-SSP-49FA**

<b>Product Name:</b>	ABT-HC-SSP-49FA
<b>Features:</b>	<ul style="list-style-type: none"> <li>• 4 swivel casters (5 3/4")</li> <li>• 115V/60 Hz</li> <li>• ETL listed</li> <li>• Exterior dimensions: 54" W x 35 3/4" D x 83 1/4" H with casters; 77 1/2" H without casters</li> <li>• 115V, 60 Hz, 9 Amps, 3/4 HP</li> <li>• ETL listed</li> <li>• Shipping weight: 495 lbs.</li> </ul>