

SPEC SHEET ABT-HC-UCFS-0120

<p>Product Name:</p>	<p>ABT-HC-UCFS-0120</p>
<p>Features:</p>	<ul style="list-style-type: none"> • Premier Countertop Freezer Freestanding • 1.7 Cu. Ft. Capacity • One solid swing door, right hinged* • Microprocessor temperature controller • Digital temperature display • Audible and visual high and low temperature alarms • Remote alarm contacts • 2 probes; 1 in air and 1 in glycol bottle • Adjustable operating temperature range: -15°C – -25°C • 2 adjustable shelves with guard rail to prevent product from being pushed against back wall • Bottom footed shelf to allow air circulation • Keyed door lock • Magnetic door gasket for positive seal • Leveling legs • Integrated door handle • Probe access port (3/8") • Powder coat exterior • High density urethane foam cabinet and door insulation • HFC-free refrigerant (R600a), foam insulation and packaging • Two year parts and labor warranty, plus an additional three year compressor parts only warranty • Pyxis®, Omnicell® and AcuDose RX® compatible • Exterior dimensions: 23 ¾" W x 24" D x 20 3/8" H • Note: This unit must have 4" clearance on sides and back for adequate ventilation • Manual defrost

SPEC SHEET **ABT-HC-UCFS-0120**

Product Name:	ABT-HC-UCFS-0120
Features:	<ul style="list-style-type: none">• 115V, 60 Hz, 1 Amp, 1/8 HP• UL/C-UL listed• ENERGY STAR CERTIFIED• Shipping weight: 106 lbs.• *This product is available left hinged, product number: ABT-HC-UCFS-0120-LH