

SPEC SHEET ABT-RFC-30M

Product Name:	ABT-RFC-30M
Features:	<ul style="list-style-type: none"> • Standard Refrigerator & Freezer Combination • 30 cu. ft. capacity (25 ref./5 freezer) • One exterior solid swing door and one interior solid swing door, right hinged • Mechanical thermostat • Refrigerator temperature range: 0°C to 10°C • Freezer temperature range -15°C to -20°C • 3 adjustable shelves, refrigerator • Keyed door lock • Magnetic door gasket for positive seal • Stainless steel door handle • Leveling legs • Two probe access ports (3/8"); one, refrigerator/one, freezer • Powder coat paint exterior • High density urethane foam cabinet and door insulation • CFC-free refrigerant (R404a), foam insulation and packaging • Two year parts and labor warranty, plus an additional three year compressor parts warranty • Exterior dimensions: 35 3/4" W x 35 1/2" D x 72 1/2" H • Manual defrost • 115V, 60 Hz, 6 Amps, 1/3 HP • Shipping weight: 383 lbs.