

**SPEC SHEET PH-ABT-HC-RFC9SS**

<b>Product Name:</b>	PH-ABT-HC-RFC9SS
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Premier Pharmacy/Vaccine Refrigerator &amp; Freezer Combination</li> <li>• 9 Cu. Ft. Capacity (4.5 Ref./4.2 Freezer)</li> <li>• Two stainless steel exterior swing doors, right hinged*</li> <li>• Microprocessor temperature controller, refrigerator and freezer</li> <li>• Digital temperature display, refrigerator and freezer</li> <li>• Adjustable operating temperature range: Refrigerator: 36°F – 46°F (2°C -8°C); Freezer: 5°F – -13°F (-15°C - -25°C)</li> <li>• Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years certification of calibration, “buffered” probe in the product simulated solution, min/max memory. F/C switchable, field installable, and visual &amp; audible temp alarm</li> <li>• Audible and visual high and low temperature alarms, refrigerator and freezer</li> <li>• Remote alarm contacts, refrigerator and freezer</li> <li>• 2 probes; 1 in air and 1 in sample bottle, in both Refrigerator and Freezer</li> <li>• 2 adjustable/1 fixed shelves refrigerator/2 fixed shelves freezer with guard rails on back in Refrigerator to prevent product from being pushed against back wall</li> <li>• Back wall evaporator cover, Refrigerator</li> <li>• Forced draft circulation, Refrigerator</li> <li>• Vaccine storage power cord warning label</li> <li>• Keyed door locks</li> <li>• Magnetic door gasket for positive seal</li> <li>• Recessed door handles</li> <li>• Leveling legs</li> <li>• Stacking kit</li> <li>• LED interior lights</li> </ul>

**SPEC SHEET PH-ABT-HC-RFC9SS**

<b>Product Name:</b>	PH-ABT-HC-RFC9SS
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Probe access port (3/8"), 2 ea. (1 for refrigerator/1 for freezer)</li> <li>• Stainless steel exterior</li> <li>• HFC-free refrigerant (R600a), foam insulation and packaging</li> <li>• Pharmacy refrigerator/freezer toolkit and temperature logs</li> <li>• Two year parts and labor warranty, plus an additional three year compressor parts warranty</li> <li>• Exterior dimensions: 23 ¾" W x 24 ½" D x 66 ¾" H</li> <li>• Manual defrost freezer / cycle defrost refrigerator</li> <li>• 115V, 60 Hz, 1.74 Amps Refrigerator, 1.5 Amps Freezer, 1/8 HP Refrigerator and Freezer</li> <li>• UL/C-UL listed</li> <li>• Shipping weight: 266 lbs.</li> <li>• *Please note this unit requires 2 outlets.</li> <li>• *This product is available left hinged, product number: PH-ABT-HC-RFC9SS-LH</li> </ul>