

SPEC SHEET PH-ABT-HC-S47G

Product Name:	PH-ABT-HC-S47G
Features:	<ul style="list-style-type: none"> • Standard Pharmacy/Vaccine Sliding Glass Door Refrigerator • 47 Cu. Ft. Capacity • Two sliding glass doors, self-closing (energy efficient, double pane vacuum insulated) • Microprocessor temperature controller • Digital temperature display • Audible and visual high and low temperature alarms • Remote alarm contacts • 2 probes; 1 in air and 1 in sample bottle • Adjustable operating temperature range: 36°F – 46°F (2°C – 8°C) • Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration, with min/max memory, °F/°C switchable • Eight adjustable shelves (adjustable in ½" increments) • Vaccine storage power cord warning label • Forced draft air circulation • Keyed door lock • Magnetic door gasket for positive seal • Integrated door handles • Swivel Casters (3 ½") • LED interior lights are safety shielded and switch controlled • Probe access port (¾") • Powder coated steel white interior and exterior • High density urethane foam cabinet insulation • Hydrocarbon, natural refrigerant (R290) • Pharmacy refrigerator/freezer toolkit and temperature logs

SPEC SHEET PH-ABT-HC-S47G

Product Name:	PH-ABT-HC-S47G
Features:	<ul style="list-style-type: none">• One year parts and labor warranty, plus an additional four year compressor parts warranty• Exterior dimensions: 54 1/8" W x 33 1/4" D x 81 3/4" H with casters; 78 1/4" H without casters• Cycle defrost• 115V, 60 Hz, 4.5 Amps. 1/4 HP• ETL listed• Shipping weight: 508 lbs.