

<b>Product Name:</b>	PH-ABT-HC-UCFS-0120A
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Premier Pharmacy/Vaccine Countertop Auto Defrost Freezer Freestanding</li> <li>• 1.3 Cu. Ft. Capacity</li> <li>• One solid swing door, right hinged*</li> <li>• Microprocessor temperature controller</li> <li>• Digital temperature display</li> <li>• Audible and visual high and low temperature alarms</li> <li>• Remote alarm contacts</li> <li>• 2 probes; 1 in air and 1 in sample bottle</li> <li>• Adjustable operating temperature range: 5°F – -13°F (-15°C – - 25°C)</li> <li>• Auto defrost cycle every 6 hours, max air temp in empty freezer: +4.34°C; max air temp in loaded freezer: -2.89°C. Max sample temp in empty freezer: -4.76°C; max sample temp in loaded freezer: - 9.58°C.</li> <li>• Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years certification of calibration, “buffered” probe in the product simulated solution, min/max memory, F/C switchable, field installable, and visual &amp; audible temp alarm</li> <li>• 2 adjustable shelves with guard rail to prevent product from being pushed against back wall</li> <li>• Bottom footed shelf to allow air circulation</li> <li>• Vaccine storage power cord warning label</li> <li>• Keyed door lock</li> <li>• Magnetic door gasket for positive seal</li> <li>• Leveling legs</li> <li>• Probe access port (3/8")</li> <li>• Powder coat exterior</li> </ul>

**SPEC SHEET PH-ABT-HC-UCFS-0120A**

<b>Product Name:</b>	PH-ABT-HC-UCFS-0120A
<b>Features:</b>	<ul style="list-style-type: none"> <li>• High density urethane foam cabinet and door insulation</li> <li>• HFC-free refrigerant (R600a), foam insulation and packaging</li> <li>• Pharmacy refrigerator/freezer toolkit and temperature logs</li> <li>• Two year parts and labor warranty, plus an additional three year compressor parts only warranty</li> <li>• Pyxis®, Omnicell® and AcuDose RX® compatible</li> <li>• Exterior dimensions: 23 ¾" W x 24" D x 20 3/8" H</li> <li>• Note: This unit must have 4" clearance on sides and back for adequate ventilation</li> <li>• Auto defrost</li> <li>• 115V, 60 Hz, 1 Amp, 1/8 HP</li> <li>• Energy Consumption 1.4 kWh/day</li> <li>• ADA compliant</li> <li>• UL/C-UL listed</li> <li>• ENERGY STAR CERTIFIED</li> <li>• Shipping weight: 106 lbs.</li> </ul>