

SPEC OF AT250B-9B


STP-JX-Tip-AT250B-9B-A2

GEB Robotic Tips for Beckman FX/NX®, Multimek AP96® and Biomek 3000®.
Product details at-a-glance

GEB Cat#	Max Volume	Tip Color	Low Retention	Sterilized	Filtered (4.3*3.5)mm	Packaging Spec.
AT250B-9B-N	250ul	Natural				96tips/rack, 10racks/pack
AT250B-9B-N-LB	250ul	Natural	Low Retention			
AT250B-9B-NS	250ul	Natural		Sterilized		5packs/case.
AT250B-9B-NS-LB	250ul	Natural	Low Retention	Sterilized		4800 tips/case,
AF250B-9B-NS	250ul	Natural		Sterilized	Filtered	50racks/case
AF250B-9B-NS-LB	250ul	Natural	Low Retention	Sterilized	Filtered	
AT250B-9B-X	250ul	Black				96tips/rack, 10racks/pack
AT250B-9B-X-LB	250ul	Black	Low Retention			
AT250B-9B-XS	250ul	Black		Sterilized		4800 tips/case,
AT250B-9B-XS-LB	250ul	Black	Low Retention	Sterilized		50racks/case
AF250B-9B-XS	250ul	Black		Sterilized	Filtered	(Conductive
AF250B-9B-XS-LB	250ul	Black	Low Retention	Sterilized	Filtered	black)
Case dimension L*W*H: 500*285*405mm.			Gross Weight : 10.2kg.			

- ✧ 250µL Liquid Level Sensing Tips
- ✧ 96 Tip/Hanging Tip Rack format
- ✧ Polypropylene
- ✧ Filtered (ith aerosol-resistant filters) /NonFiltered
- ✧ Presterile / Nonsterile
- ✧ Certified RNase-/DNase-free
- ✧ Nonpyrogenic
- ✧ Endotoxin-Safe meets USP/FDA standards
- ✧ Comply With USP 3127 (Heavy Metal Contamination) standards
- ✧ Compatibility: Beckman FX/NX®, Multimek AP96® and Biomek 3000® workstations.

Drawing: (No tolerance Diameter tolerance : $\pm 0.05\text{mm}$, Height tolerance : $\pm 0.25\text{mm}$)



Prepared by: Anne Mar.26, 2021 Approve: WeiMing Yu Mar.26,2021