

Diamond RV-Pette-PRO™

COMPATIBILITY CHART

Dispenser Syringe Tips

Wide range of sizes match all sizes offered by other manufacturers and eliminates the need to change protocols or practice!

12
volume options!

Diamond RV-Pette-PRO™ Dispenser Syringe Tip Size:

Pipettor Brand/Model:	0.05mL	0.1mL	0.2mL	0.5mL	1.0mL	1.25mL	2.5mL	5.0mL	10.0mL	12.5mL	25mL	50mL
Diamond RV-Pette™	YES	YES	YES	YES	YES	YES	YES	YES	YES	YES	YES	YES
Diamond RV-Pette PRO™	YES	YES	YES	YES	YES	YES	YES	YES	YES	YES	YES	YES
Bel-Art Roxy-M™ Repeating Pipettor	✓	—	—	✓	—	✓	—	✓	—	✓	✓	✓
Capp Rhythm™ Repeater	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Capp R10™	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Brand® Handystep® S	—	✓	—	✓	✓	✓	✓	✓	✓	✓	✓	✓
Brand® Handystep® Electronic	—	✓*	—	✓*	✓*	✓*	✓*	✓*	✓*	✓*	✓*	✓*
Biohit® Mechanical Stepper	✓	—	—	✓	—	✓	✓	✓	—	✓	✓	✓
eLine® Light Dispenser	—	✓	✓	✓	✓	—	✓	✓	✓	—	✓	✓
eLine® Pro Dispenser	—	✓	✓	✓	✓	—	✓	✓	✓	—	✓	✓
Eppendorf® Repeater 4780	✓	—	—	✓	—	✓	✓	✓	—	✓	✓	✓
Fisherbrand™	—	✓	—	✓	✓	✓	✓	✓	✓	✓	✓	✓
Gilson® DISTRIMAN®	—	—	—	—	—	✓	—	—	—	✓	—	—
Gilson® Repetman®	—	✓*	—	✓*	✓*	✓*	✓*	✓*	✓*	✓*	✓*	✓*
Nichiryo NichiMate	✓	—	—	✓	—	✓	✓	✓	—	✓	✓	✓
Rainin® AutoRep S	—	✓	—	✓	✓	✓	✓	✓	✓	✓	✓	✓
Rainin® AutoRep E	—	✓*	—	✓*	✓*	✓*	✓*	✓*	✓*	✓*	✓*	✓*
VWR® Repeating Dispenser	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

*Tip-size needs to be entered manually.

MIDSCI™
Laboratory Equipment & Supplies

All instruments and names of manufacturers mentioned herein are the trademarks of their respective companies and are used for reference purposes only. MIDSCI expressly disclaims any affiliation to them.

ORDER TODAY! midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998