

SPEC OF AT1000R-B/9/E

STP-JX-Tip-AT1000R-B/9/E-A2

GEB Robotic Tips for Rainin® can be used in proteomics, drug development, genomics and other fields.

Product details at-a-glance

GEB Cat#	Max Volume	Tip Color	Low Retention	Sterilized	Filtered (6.6*3.5)mm	Packaging Spec.
AT01MR-9-N	1,000ul	Natural				96 tips/rack, 4,800 tips/case
AT01MR-9-NS	1,000ul	Natural		Sterilized		
AF01MR-9-NS	1,000ul	Natural		Sterilized	Filtered	
AF01MR-9-NS-LB	1,000ul	Natural	Low Retention	Sterilized	Filtered	
Case dimension L*W*H: 595*285*475mm			Gross Weight : 13.5kg			

- ✧ 1000µL Natural Tips
- ✧ 96 Tip/Hanging Tip Rack format/Blut
- ✧ Polypropylene
- ✧ Certified DNase/RNase free
- ✧ Nonpyrogenic
- ✧ Presterile/Nonsterile
- ✧ Filtered(Filtered With aerosol-resistant filters)/NonFiltered
- ✧ Endotoxin-Safe meets USP/FDA standards
- ✧ Comply With USP 3127 (Heavy Metal Contamination) standards
- ✧ Compatibility: Rainin® liquid workstations.

Drawing: (No tolerance Diameter tolerance : ±0.05mm , Height tolerance : ±0.25mm)

