

SPEC OF AT200R-B/9/E

STP-JX-Tip-AT200R-B/9/E-A2

GEB Robotic Tips for Rainin® can be used in proteomics, drug development, genomics and other fields.

Product details at-a-glance

GEB Cat#	Max Volume	Tip Color	Low Retention	Sterilized	Filtered (4.5*3.5)mm	Packaging Spec.
AT200R-9-N	200ul	Natural				96 tips/rack, 10 racks/pack, 5 packs/case.
AT200R-9-NS	200ul	Natural		Sterilized		
AF200R-9-NS	200ul	Natural		Sterilized	Filtered	50 racks/case, 4,800 tips/case.
AF200R-9-NS-LB	200ul	Natural	Low Retention	Sterilized	Filtered	
Case dimension L*W*H: 500*285*405mm			Gross Weight : 10.0 kg			

- ✧ 200µL Natural Tips
- ✧ 96 Tip/Hanging Tip Rack format/Blut
- ✧ Polypropylene
- ✧ Certified DNase/RNase free
- ✧ Nonpyrogenic
- ✧ Presterile/Nonsterile
- ✧ Filtered(Filtered With aerosol-resistant filters)/NonFiltered
- ✧ Endotoxin-Safe meets USP/FDA standards
- ✧ Comply With USP 3127 (Heavy Metal Contamination) standards
- ✧ Compatibility: Rainin® liquid workstations

Drawing: (No tolerance Diameter tolerance : $\pm 0.05\text{mm}$, Height tolerance : $\pm 0.25\text{mm}$)

