

COMFORT AND DEXTERITY FOR EVERYDAY HAND PROTECTION.



TF-065-095-NT
POWDER-FREE LATEX EXAMINATION GLOVES


PROVIDES BARRIER PROTECTION AGAINST CHEMICALS, PATHOGENS AND MORE FOR MAXIMIZE DEXTERITY

- Natural latex provides advanced chemical protection
- Exceeds 1.5 AQL laboratory standards





Optimal thickness for durability and tactile sensitivity

Natural latex provides advanced chemical protection

Tear-resistant beaded cuff for easy donning



TF TF-065-095-NT

WEIGHT	LENGTH	THICKNESS	COLOR
 6.5g (large size)	 9.5"	 4.7 mil	 NATURAL





TF-065-095-NT POWDER-FREE LATEX EXAMINATION GLOVES



PRODUCT DESCRIPTION

SPECIFICATIONS

PRODUCT	TrueForm
GRADE	Powder-Free Exam
MATERIAL	Latex
THICKNESS FINGER	5.5mil 0.14mm
PALM	4.7mil 0.12mm
COLOR	Natural
SURFACE	Fully Textured
LENGTH	9.5" 242mm
WEIGHT	6.5g
FIT TYPE	Ambidextrous
CUFF TYPE	Beaded
WATERTIGHT AQL	1.5
ORIGIN COUNTRY	Malaysia

PHYSICAL PROPERTIES

	PRE-AGING	POST-AGING
TENSILE STRENGTH:	19 MPA	16 MPA
ELONGATION:	850%	750%



* Specifications based on size Large glove

STANDARDS & CERTIFICATIONS

EN 420-2003
EN455-2000-2015
ASTM D3578
ASTM F1671-97

FDA 21 CFR 177.2600 Component materials comply with FDA regulations for food contact
SHA Dermatologically Approved
FDA 510K Clearance

CE 0493

EN ISO 374-1:2016/TYPE C EN ISO 374-5:2016



Please visit swsafety.com/symbols for information on standards used on this document.

ORDERING

ITEM NO.	SIZE	PALM WIDTH	PACKAGING
L006201	XS 5½ – 6	3.0" 74mm	100 gloves per box 10 boxes per case
L006202	S 6½ – 7	3.5" 87mm	100 gloves per box 10 boxes per case
L006203	M 7½ – 8	3.7" 95mm	100 gloves per box 10 boxes per case
L006204	L 8½ – 9	4.3" 108mm	100 gloves per box 10 boxes per case

APPLICATIONS

- Chemical Handling
- Dairy
- Healthcare
- Industrial
- Janitorial Maintenance
- Laboratory
- Life Sciences
- Tattooing and Cosmetics

Keep out of sunlight. Store in a cool, dry place.
Keep away from sources of ozone or ignition.



R_L00620X0522

