

STRIPSPIN™ 12D

Digital Mini Centrifuge

MIDSCI™
Laboratory Equipment & Supplies



- Unique rotor accepts 12-position PCR strips
- Digital control up to 5,500 rpm
- “Mini” sized, less than 8 inches wide
- Timer up to 99 mins.



Holds up to 4 x PCR strips or 48 x 0.2ml tubes

The StripSpin 12D from Benchmark is the only mini centrifuge designed for the spin down of 12-position PCR strips. The centrifuge features digital speed control up to 5,500 rpm, time control up to 99 min., and a safety interlock to stop the centrifuge if the lid is opened during operation.

With built in rcf conversion, speed can be set/displayed in rpm or rcf. The StripSpin 12D is the perfect solution for centrifugation of both 8 and 12-position PCR strips (or 48 x 0.2ml tubes) and includes a 2 year warranty.

Technical Data

Speed:	5,500rpm
Max. g-force:	1,690 xg
Capacity:	4 x 12 position PCR strips 4 x 8 position PCR strips 48 x 0.2ml PCR tubes
Dimensions: (W x D x H)	7.6 x 9 x 4.7 in 19.4 x 23 x 12 cm
Weight:	4.4 lbs / 2 kg
Electrical:	100 to 240V, 50-60Hz
Warranty:	2 Years

Also available...



Item: C2212



Includes COMBI-Rotor™
for tubes and strips

Ordering Information

C2248* **StripSpin** Mini Centrifuge, with rotor for up to 4 x PCR strips (8 or 12 position) or 48 x 0.2ml tubes

*Includes US plug. For EU plug, add (-E) to the item number.

Benchmark
Scientific

ORDER TODAY! midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998