

## flip&grip pipette holder

### Flexible bench or under-shelf storage

Unique pipette holder that gives you the option of using on the bench top, or hanging under shelf holder to help maximize your working space.

- Will hold up to 12 single channel pipette, or a combo of 3 single and 3 multichannel pipettes
- Holds most common brands
- Rubber grips prevent scratching pipette
- Holder rotates for easy pipette selection and storage
- Weighted base
- Fixing screw provided for under shelf option



Pipettes can be held in the rubber grips and/or the plunger hooks



Mount base under shelf or counter

### Specifications

H x W Bench Top Style	H x W Under Shelf Style	Color
332 x 193 mm	187 x 193 mm	Blue ■

### Ordering Information

ITEM #	DESCRIPTION
HS120026	Flip and Grip Pipette Stand, holds 12 single channel or 3 multichannel and 3 single

Rack rotates for easy pipette selection