

July 14, 2022

To whom it may concern:

Labcon is conducting a recall of the following products:

7626-500-000-9

and

7626-540-000-9

Product description: U-pette dispenser

Lot numbers affected: any lot number with 2 in the position as highlighted below.

XXXC8-2XX

Purchase orders that are possibly affected: please see the list below

119387

Reason for the recall: While performing a product audit it was discovered that metal cannulas in the U-pette dispensers could occasionally get displaced when the dispenser is inserted in the Vacutainer stopper.

Action required: please check your inventory for the presence of the lot numbers mentioned above. Please contact the end-users if needed.

Supplier contact for credit and/or replacement: If any of these lots are found in your stock please contact custsupport@labcon.com

CAPA statement: A new testing protocol is implemented for U-pette cannula testing as of 07/14/2022

Please contact custsupport@labcon.com if you have any questions or concerns.

Lily Remennik
Customer Support/QA Manager