AbsorbMax™

For use with black microplates



Application:

Fluorescence assays, light sensitive sample protection



Temperature Range:

-40°C to +80°C

CATALOG#	DESCRIPTION	QUANTITY
BK-50	AbsorbMax	50 Sheets



Features and Benefits

- Black vinyl film that can be applied to either the top or bottom of the plate
- Minimizes background fluorescence and cross talk between plate wells
- Ideal as a back seal for black plates with clear well bottoms for top reading measurements
- Absorbs stray light, reducing the background caused by light scatter
- Protects light sensitive samples from degradation
- Features a single split liner for easy application

Bundling Opportunities

- Black microplates
- Pipettors
- Pipet tips
- Reservoirs



