# FoilSeal™ RollMate™

Pre-scored foil in an easy-to-use roll format



## **Application:**


PCR, light sensitive assays, short-term storage and robotics



### **Temperature Range:**

-20°C to +120°C

CATALOG#	DESCRIPTION	QUANTITY
F2-RL-500	FoilSeal RollMate	500/Roll



#### **Features and Benefits**

- · Certified DNase- & RNase- free
- Pre-scored film on a roll is a more economical and efficient alternative to rapidly seal microplates; faster than applying single sheeted film
- Self-standing, non-slip supports hold the pre-scored film roll securely in place for quick and easy benchtop dispensing
- The roll unwind tension can be easily adjusted to user preference
- Continuous backing liner can be torn off or fed directly into a wastebasket for disposal
- Pierceable with a pipet tip or robotic probe for direct sample recovery
- 2.0mil aluminum foil removes easily without tearing
- The conformable nature of foil allows for the visual verification of the seal's security around each well
- Not recommended for long-term cold storage

## **Bundling Opportunities**

- PCR plates
- PCR reagents and buffers
- Thermal cyclers
- Pipettors and pipet tips
- Reservoirs



