

**SPEC SHEET ABT-HC-16S**

<b>Product Name:</b>	ABT-HC-16S
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Premier Solid Door Laboratory Refrigerator</li> <li>• 16 Cu. Ft. capacity</li> <li>• One swing solid door, self-closing, right hinged, not reversible</li> <li>• Microprocessor Temperature Controller</li> <li>• Adjustable operating temperature range: 1°-10°C</li> <li>• Standard Digital Temperature Display and Alarm Module with battery back-up to support temperature monitoring and alarms</li> <li>• Sample simulated digital temperature display; .1° Resolution</li> <li>• Adjustable temperature display in Celsius or Fahrenheit</li> <li>• Audible and visual high/low temperature alarms</li> <li>• 2 probes; 1 in air and 1 in sample bottle</li> <li>• Sensor error alarm, Power failure alarm</li> <li>• Resettable 1-point Min/Max History</li> <li>• Remote alarm contacts</li> <li>• Alarm validation (audible, visual, and remote alarm)</li> <li>• Alarm mute with ring back</li> <li>• Three adjustable shelves (adjustable in 1/2" increments)</li> <li>• Forced draft air circulation</li> <li>• Keyed door lock</li> <li>• Magnetic door gasket for positive seal</li> <li>• Heavy duty door handle</li> <li>• LED interior lights are safety shielded and switch controlled</li> <li>• Swivel casters (3 1/2") probe access port (3/8")</li> <li>• White powder-coated steel exterior and interior</li> <li>• Hydrocarbon, natural refrigerant (R290)</li> <li>• Two year parts and labor warranty, plus an additional three compressor parts warranty</li> <li>• Exterior dimensions: 25" W x 29 3/4" D x 79" H with casters; 75 1/2" H without casters</li> <li>• Cycle defrost</li> <li>• 115V, 60 Hz, 3 Amps, 1/5 HP</li> <li>• ETL listed</li> <li>• Shipping weight: 277 lbs.</li> </ul>