

SPEC SHEET **ABT-HC-36S**

| | |
|----------------------|--|
| Product Name: | ABT-HC-36S |
| Features: | <ul style="list-style-type: none"> • Premier Solid Door Laboratory Refrigerator • 36 Cu. Ft. capacity • Two swing solid doors, self-closing, right hinged, not reversible • Microprocessor Temperature Controller • Adjustable operating temperature range: 1°-10°C • Standard Digital Temperature Display and Alarm Module with battery back-up to support temperature monitoring and alarms • Sample simulated digital temperature display; .1° Resolution • Adjustable temperature display in Celsius or Fahrenheit • Audible and visual high/low temperature alarms • 2 probes; 1 in air and 1 in sample bottle • Sensor error alarm, Power failure alarm • Resettable 1-point Min/Max History • Remote alarm contacts • Alarm validation (audible, visual, and remote alarm) • Alarm mute with ring back • Eight adjustable shelves (adjustable in ½" increments) • Forced draft air circulation • Keyed door lock • Magnetic door gasket for positive seal • Heavy duty door handle • LED interior lights are safety shielded and switch controlled |

SPEC SHEET **ABT-HC-36S**

| | |
|-----------------------------|---|
| <p>Product Name:</p> | <p>ABT-HC-36S</p> |
| <p>Features:</p> | <ul style="list-style-type: none"> • Swivel casters (3 ½"), probe access port (3/8") • White powder-coated steel exterior and interior • Hydrocarbon, natural refrigerant (R290) • Two year parts and labor warranty, plus an additional three compressor parts warranty • Exterior dimensions: 39 5/8" W x 34 3/4" D x 81 3/4" H with casters; 78 1/4" H without casters • Cycle defrost • 115V, 60 Hz, 4.5 Amps, 1/4 HP • ETL listed • Shipping weight: 360 lbs. |