

SPEC SHEET **ABT-HC-CP-47**

Product Name:	ABT-HC-CP-47
Features:	<ul style="list-style-type: none"> • Premier Glass Door Chromatography Refrigerator • 47 Cu. Ft. capacity • Two sliding glass doors, self-closing (energy efficient, double pane vacuum insulated) • Microprocessor Temperature Controller • Adjustable operating temperature range: 1°-10°C • Digital Temperature Display and Alarm Module • Sample simulated digital temperature display; .1° Resolution • 2 probes; 1 in air and 1 in sample bottle • Adjustable temperature display in Celsius or Fahrenheit • Audible and visual high/low temperature alarms • Sensor error alarm; Power failure alarm • Resettable 1-point Min/Max History • Remote alarm contacts • Alarm validation (audible, visual, and remote alarm) • Alarm mute with ring back • Eight adjustable shelves (adjustable in 1/2" increments) • Forced draft air circulation • Keyed door lock • One GFCI protected duplex electrical outlet • Two 2" access ports with spring-loaded, self-closing gasketed metal covers • Magnetic door gasket for positive seal and integrated door handles • LED interior lights are safety shielded and switch controlled • Swivel Casters (3 1/2") Probe access port (3/4")

SPEC SHEET ABT-HC-CP-47

Product Name:	ABT-HC-CP-47
Features:	<ul style="list-style-type: none"> • Advanced hydrocarbon detector • White powder-coated steel exterior and interior • High density urethane foam cabinet insulation • Hydrocarbon, natural refrigerant (R290) • Two year parts and labor warranty, plus an additional three compressor parts warranty • Exterior dimensions: 55 5/8" W* x 33 1/4" D x 81 3/4" H with casters; 78 1/4" H without casters. *Width w/access ports closed;63 5/8" width w/access ports open • Cycle defrost • 115V, 60 Hz, 4.5 Amps, 1/4 HP • ETL listed • Shipping weight: 508 lbs. • Note: Unit requires two (2) outlets