

**SPEC SHEET**    **ABT-HC-CP-72**

<b>Product Name:</b>	ABT-HC-CP-72
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Premier Glass Door Chromatography Refrigerator</li> <li>• 72 Cu. Ft. capacity</li> <li>• Three swing glass doors, self-closing, right hinged, not reversible</li> <li>• Microprocessor Temperature Controller</li> <li>• Adjustable operating temperature range: 1°-10°C Digital</li> <li>• Temperature Display and Alarm Module</li> <li>• Sample simulated digital temperature display; .1° Resolution</li> <li>• Adjustable temperature display in Celsius or Fahrenheit</li> <li>• Audible and visual high/low temperature alarms</li> <li>• 2 probes; 1 in air and 1 in sample bottle</li> <li>• Sensor error alarm, Power failure alarm</li> <li>• Resettable 1-point Min/Max History</li> <li>• Remote alarm contacts</li> <li>• Alarm validation (audible, visual, and remote alarm)</li> <li>• Alarm mute with ring back</li> <li>• Twelve adjustable shelves (adjustable in ½" increments)</li> <li>• Forced draft air circulation</li> <li>• Keyed door lock</li> <li>• One GFCI protected duplex electrical outlet</li> <li>• Two 2" access ports with spring-loaded, self-closing gasketed metal covers</li> <li>• Magnetic door gasket for positive seal</li> <li>• Heavy duty door handle</li> <li>• LED interior lights are safety shielded and switch controlled</li> <li>• Swivel casters (3 ½"), Probe access port (3/4")</li> </ul>

**SPEC SHEET ABT-HC-CP-72**

<b>Product Name:</b>	ABT-HC-CP-72
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Advanced hydrocarbon detector</li> <li>• White powder-coated steel exterior and interior</li> <li>• High density urethane foam cabinet insulation</li> <li>• Hydrocarbon, natural refrigerant (R290)</li> <li>• Two year parts and labor warranty, plus an additional three year compressor parts warranty</li> <li>• Exterior dimensions: 82 1/2" W* x 34" D x 78 1/4" H without casters; 81 3/4" H with casters. *Width w/access ports closed; 90 1/2" width w/access ports open</li> <li>• Cycle defrost</li> <li>• 115V, 60 Hz, 7.5 Amps, 1/2 HP</li> <li>• ETL listed</li> <li>• Shipping weight: 743 lbs.</li> <li>• Note: This unit requires two (2) outlets</li> </ul>