

SPEC SHEET ABT-HC-CS-26

Product Name:	ABT-HC-CS-26
Features:	<ul style="list-style-type: none"> • Standard Glass Door Chromatography Refrigerator • 26 Cu. Ft. capacity • One swing glass door, self-closing, right hinged, not reversible • Microprocessor temperature controller • Digital temperature display • Audible and visual high and low temperature alarms • 2 probes; 1 in air and 1 in sample bottle • Remote alarm contacts • Adjustable operating temperature range: 1°C – 10°C • Four adjustable shelves (adjustable in ½" increments) • Forced draft air circulation • Keyed door lock • One GFCI protected duplex electrical outlet • One 2" access port with spring-loaded, self-closing gasketed metal cover • Magnetic door gasket for positive seal • Heavy duty door handle • LED interior lights are safety shielded and switch controlled • Swivel Casters (3 ½") Probe access port (3/8") • Advanced hydrocarbon detector • White powder-coated steel exterior and interior • High density urethane foam cabinet insulation • Hydrocarbon, natural refrigerant (R290) • Two year parts and labor warranty, plus an additional three year compressor parts warranty

SPEC SHEET **ABT-HC-CS-26**

Product Name:	ABT-HC-CS-26
Features:	<ul style="list-style-type: none">• Exterior dimensions: 29 1/8" W x 36 1/2" D x 78 1/4" H without casters; 81 3/4" H with casters• Cycle defrost• 115V, 60 Hz, 3 Amps, 1/5 HP• ETL listed• Shipping weight: 355 lbs.• Note: This unit requires two (2) outlets