

SPEC SHEET ABT-HC-CS-36

Product Name:	ABT-HC-CS-36
Features:	<ul style="list-style-type: none"> • Standard Glass Door Chromatography Refrigerator • 36 Cu. Ft. capacity • Two swing glass doors, self-closing, right hinged, not reversible • Microprocessor temperature controller • Digital temperature display • Audible and visual high and low temperature alarms • 2 probes; 1 in air and 1 in sample bottle • Remote alarm contacts • Adjustable operating temperature range: 1°C – 10°C • Keyed door locks • One GFCI protected duplex electrical outlet • One 2" access port with spring-loaded, self-closing gasketed metal cover • Magnetic door gasket for positive seal • Eight adjustable shelves (adjustable in ½" increments) • Forced draft air circulation • Heavy duty door handles • LED interior lights are safety shielded and switch controlled • Swivel Casters (3 ½") • Probe access port (3/8") • Advanced hydrocarbon detector • White powder-coated steel exterior and interior • High density urethane foam cabinet insulation • Hydrocarbon, natural refrigerant (R290) • Two year parts and labor warranty, plus an additional three year compressor parts warranty

SPEC SHEET ABT-HC-CS-36

Product Name:	ABT-HC-CS-36
Features:	<ul style="list-style-type: none"> • Exterior dimensions: 40 3/8" W x 34 3/4" D x 81 3/4" H with casters; 78 1/4" H without casters • Cycle defrost • 115V, 60 Hz, 4.5 Amps, 1/4 HP • ETL listed • Shipping weight: 405 lbs. • Note: This unit requires two (2) outlets