

**SPEC SHEET ABT-HC-FFP-14**

<b>Product Name:</b>	ABT-HC-FFP-14
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Standard Natural Refrigerant Flammable Storage Freezer</li> <li>• 14 Cu. Ft. Capacity</li> <li>• One solid swing door, right hinged</li> <li>• Microprocessor temperature controller</li> <li>• Digital temperature display with sample simulated probe</li> <li>• Audible and visual high and low temperature alarms</li> <li>• Remote alarm contacts</li> <li>• Adjustable operating temperature range: -15°C – -25°C</li> <li>• Six fixed shelves</li> <li>• Keyed door lock</li> <li>• Magnetic door gasket for positive seal</li> <li>• Leveling legs</li> <li>• Probe access port (3/8")</li> <li>• Powder coat exterior</li> <li>• High density urethane foam cabinet and door insulation</li> <li>• Natural, Hydrocarbon Refrigerant (R290)</li> <li>• Two year parts and labor warranty</li> <li>• Utilizes ATEX and non-sparking components for safety</li> <li>• Exterior dimensions: 23 5/8 "W x 24 3/8" D x 73"H</li> <li>• Note: This unit must have 4" clearance on sides and back for adequate ventilation</li> <li>• Manual defrost</li> <li>• 115V, 60 Hz, 3.5 Amps, 1/3 HP</li> <li>• UL/C-UL or ETL Listed</li> <li>• Shipping weight: 195 lbs.</li> </ul>