

SPEC SHEET ABT-HC-FRP-17

Product Name:	ABT-HC-FRP-17
Features:	<ul style="list-style-type: none"> • Standard Natural Refrigerant Flammable Storage Refrigerator • 17 Cu. Ft. capacity • One solid swing door, right hinged, not reversible • Microprocessor Temperature Controller • Adjustable operating temperature range: 1°-10°C • Digital Temperature Display and Alarm Module with battery back-up to support temperature monitoring and alarms • Sample simulated digital temperature display; .1° Resolution • Adjustable temperature display in Celsius or Fahrenheit • Audible and visual high/low temperature alarms • Sensor error alarm • Power failure alarm • Resettable 1-point Min/Max History • Remote alarm contacts • Alarm validation (audible, visual, and remote alarm) • Alarm mute with ring back • 5 adjustable shelves • Keyed door lock • Magnetic door gasket for positive seal • Heavy duty door handle • Leveling legs • One probe access port (3/8") • White textured cabinet and door • High density urethane foam cabinet and door insulation • Natural, Hydrocarbon Refrigerant (R600a)

SPEC SHEET **ABT-HC-FRP-17**

Product Name:	ABT-HC-FRP-17
Features:	<ul style="list-style-type: none">• Two year parts and labor warranty• Exterior dimensions: 30 3/4" W x 28 5/8"D x 67 3/4"H• Cycle defrost• 115V, 60 Hz, 3 Amps, 1/4 HP• ETL listed• Units meet NFPA guidelines 45 and 70• Shipping weight: 240 lbs.