

**SPEC SHEET ABT-HC-LP-33**

<b>Product Name:</b>	ABT-HC-LP-33
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Premier Glass Door Laboratory Refrigerator</li> <li>• 33 Cu. Ft. capacity</li> <li>• Two sliding glass doors, self-closing (energy efficient, double pane vacuum insulated)</li> <li>• Microprocessor Temperature Controller</li> <li>• Adjustable operating temperature range: 1°-10°C</li> <li>• Digital Temperature Display and Alarm Module</li> <li>• Sample simulated digital temperature display; .1° Resolution</li> <li>• Adjustable temperature display in Celsius or Fahrenheit</li> <li>• Audible and visual high/low temperature alarms</li> <li>• 2 probes; 1 in air and 1 in sample bottle</li> <li>• Sensor error alarm; Power failure alarm</li> <li>• Resettable 1-point Min/Max History</li> <li>• Remote alarm contacts</li> <li>• Alarm validation (audible, visual, and remote alarm)</li> <li>• Alarm mute with ring back</li> <li>• Eight adjustable shelves (adjustable in ½" increments)</li> <li>• Keyed door lock</li> <li>• Forced draft air circulation</li> <li>• Magnetic door gasket for positive seal</li> <li>• Integrated door handles</li> <li>• LED interior lights are safety shielded and switch controlled</li> <li>• Swivel Casters (3 ½"), Probe access port (3/8")</li> <li>• White powder-coated steel exterior and interior</li> <li>• High density urethane foam cabinet insulation</li> </ul>

**SPEC SHEET**    **ABT-HC-LP-33**

<b>Product Name:</b>	ABT-HC-LP-33
<b>Features:</b>	<ul style="list-style-type: none"><li>• Hydrocarbon, natural refrigerant (R290)</li><li>• Two year parts and labor warranty, plus an additional three compressor parts warranty</li><li>• Exterior dimensions: 39 5/8" W x 33 1/4" D x 81 3/4" H with casters; 78 1/4" H without casters</li><li>• Cycle defrost</li><li>• 115V, 60 Hz, 4.5 Amps, 1/4 HP</li><li>• ETL listed</li><li>• Shipping weight: 419 lbs.</li></ul>