

SPEC SHEET ABT-HC-LP-47

Product Name:	ABT-HC-LP-47
Features:	<ul style="list-style-type: none"> • Premier Glass Door Laboratory Refrigerator • 47 Cu. Ft. capacity • Two sliding glass doors, self-closing (energy efficient, double pane vacuum insulated) • Microprocessor Temperature Controller • Adjustable operating temperature range: 1°-10°C • Digital Temperature Display and Alarm Module • Sample simulated digital temperature display; .1° Resolution • 2 probes; 1 in air and 1 in sample bottle • Adjustable temperature display in Celsius or Fahrenheit • Audible and visual high/low temperature alarms • Sensor error alarm; power failure alarm • Resettable 1-point Min/Max History • Remote alarm contacts • Alarm validation (audible, visual, and remote alarm) • Alarm mute with ring back • Eight adjustable shelves (adjustable in ½" increments) • Forced draft air circulation • Keyed door lock • Magnetic door gasket for positive seal • Integrated door handles • LED interior lights are safety shielded and switch controlled • Swivel Casters (3 ½"); Probe access port (⅜") • White powder-coated steel exterior and interior • High density urethane foam cabinet insulation • Hydrocarbon, natural refrigerant (R290)

SPEC SHEET **ABT-HC-LP-47**

Product Name:	ABT-HC-LP-47
Features:	<ul style="list-style-type: none">• Two year parts and labor warranty, plus an additional three compressor parts warranty• Exterior dimensions: 54 1/8" W x 33 1/4" D x 81 3/4" H with casters; 78 1/4" H without casters• Cycle defrost• 115V, 60 Hz, 4.5 Amps, 1/4 HP• ETL listed• Shipping weight: 508 lbs.