

SPEC SHEET ABT-HC-MFP-14

Product Name:	ABT-HC-MFP-14
Features:	<ul style="list-style-type: none"> • Standard Natural Refrigerant Manual Defrost Laboratory Freezer • 14 Cu. Ft. Capacity • One solid swing door, right hinged • Microprocessor temperature controller and Digital temperature display • Audible and visual high and low temperature alarms • Remote alarm contacts • 2 probes; 1 in air and 1 in glycol bottle • Adjustable operating temperature range: -15°C – -25°C • Six fixed shelves • Keyed door lock and Magnetic door gasket for positive seal • Leveling legs • Probe access port (3/8") • Powder coat exterior • High density urethane foam cabinet and door insulation • Natural, Hydrocarbon Refrigerant (R290) • Two year parts and labor warranty, plus additional three year compressor parts warranty • Exterior dimensions: 23 5/8 "W x 24 3/8" D x 73" H • Note: This unit must have 4" clearance on sides and back for adequate ventilation • Manual defrost • 115V, 60 Hz, 3.5 Amps, 1/3 HP • ENERGY STAR CERTIFIED • UL/C-UL or ETL Listed • Shipping weight: 195 lbs.