

**SPEC SHEET ABT-HC-RFP-14**

<b>Product Name:</b>	ABT-HC-RFP-14
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Standard Natural Refrigerant Laboratory Refrigerator 14 Cu. Ft. capacity</li> <li>• One solid swing door, right hinged, not reversible</li> <li>• Microprocessor Temperature Controller</li> <li>• Adjustable operating temperature range: 1°-10°C</li> <li>• Digital Temperature Display and Alarm Module with battery back-up to support temperature monitoring and alarms</li> <li>• Sample simulated digital temperature display; .1° Resolution</li> <li>• Adjustable temperature display in Celsius or Fahrenheit</li> <li>• Audible and visual high/low temperature alarms</li> <li>• Sensor error alarm</li> <li>• Power failure alarm</li> <li>• Resettable 1-point Min/Max History</li> <li>• Remote alarm contacts</li> <li>• Alarm validation (audible, visual, and remote alarm)</li> <li>• Alarm mute with ring back</li> <li>• 5 adjustable shelves</li> <li>• Keyed door lock</li> <li>• Magnetic door gasket for positive seal</li> <li>• Heavy duty door handle</li> <li>• Leveling legs</li> <li>• One probe access port (3/8")</li> <li>• White textured cabinet and door</li> <li>• High density urethane foam cabinet and door insulation</li> <li>• Natural, Hydrocarbon Refrigerant (R600a)</li> <li>• Two year parts and labor warranty</li> </ul>

**SPEC SHEET**    **ABT-HC-RFP-14**

<b>Product Name:</b>	ABT-HC-RFP-14
<b>Features:</b>	<ul style="list-style-type: none"><li>• Exterior dimensions: 23 5/8" W x 24 3/8" D x 73" H</li><li>• Cycle defrost</li><li>• 115V, 60 Hz, 3 Amps, 1/5 HP</li><li>• ETL listed</li><li>• Shipping weight: 200 lbs.</li></ul>