

SPEC SHEET ABT-HC-UCBI-0404G

Product Name:	ABT-HC-UCBI-0404G
Features:	<ul style="list-style-type: none"> • Premier Undercounter Refrigerator Built-In • 4.6 Cu. Ft. Capacity • One glass swing door, right hinged* • Microprocessor temperature controller • Digital temperature display • Audible and visual high and low temperature alarms • Remote alarm contacts • 2 probes; 1 in air and 1 in glycol bottle • Adjustable operating temperature range: 1°C – 10°C • Three shelves (two adjustable/one fixed) with guard rail on back to prevent product from pushing against back wall • Back wall evaporator cover • Forced draft circulation • Keyed door lock • Magnetic door gasket for positive seal • Leveling legs • LED interior light • Probe access port (3/8") • Powder coat exterior • High density urethane foam cabinet and door insulation • HFC-free refrigerant (R600a), foam insulation and packaging • Two year parts and labor warranty, plus an additional three year compressor parts only warranty • Pyxis®, Omnicell® and AcuDose RX® compatible • Exterior dimensions: 23 ¾" W x 25 ½" D x 33 3/8" H • Energy Star Certified for Laboratory Refrigerators and Freezers

SPEC SHEET **ABT-HC-UCBI-0404G**

Product Name:	ABT-HC-UCBI-0404G
Features:	<ul style="list-style-type: none">• Cycle defrost• 115V, 60 Hz, 1.74 Amps, 1/8 HP• Energy Consumption 1.15 kWh/day UL/C-UL listed• ENERGY STAR CERTIFIED• Shipping weight: 150 lbs.• *This product is available left hinged, product number: ABT-HC-UCBI-0404G-LH