

**SPEC SHEET ABT-HC-UCFS-0104**

<b>Product Name:</b>	ABT-HC-UCFS-0104
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Premier Countertop Freestanding Refrigerator</li> <li>• 1 Cu. Ft. Capacity</li> <li>• One solid swing door, right hinged*</li> <li>• Microprocessor temperature controller</li> <li>• Digital temperature display</li> <li>• Audible and visual high and low temperature alarms</li> <li>• Remote alarm contacts</li> <li>• 2 probes; 1 in air and 1 in glycol bottle</li> <li>• Adjustable operating temperature range: 1°C – 10°C</li> <li>• 2 adjustable shelves</li> <li>• Back wall evaporator cover</li> <li>• Forced draft circulation</li> <li>• Keyed door lock</li> <li>• Integrated door handle</li> <li>• Magnetic door gasket for positive seal</li> <li>• Leveling legs</li> <li>• LED interior light</li> <li>• Probe access port (3/8")</li> <li>• Powder coat exterior</li> <li>• High density urethane foam cabinet and door insulation</li> <li>• HFC-free refrigerant (R600a), foam insulation and packaging</li> <li>• Two year parts and labor warranty, plus an additional three year compressor parts only warranty</li> <li>• Pyxis®, Omnicell® and AcuDose RX® compatible</li> <li>• Exterior dimensions: 17 ¼" W x 19 ¼" D x 21 ¼" H</li> <li>• Note: This unit must have 4" clearance on sides and back for adequate ventilation</li> </ul>

**SPEC SHEET**    **ABT-HC-UCFS-0104**

<b>Product Name:</b>	ABT-HC-UCFS-0104
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Energy Star Certified for Laboratory Refrigerators and Freezers</li> <li>• Cycle defrost</li> <li>• 115V, 60 Hz, .90 Amps, 1/8</li> <li>• Energy Consumption 0.71 kWh/day UL/C-UL listed</li> <li>• ENERGY STAR CERTIFIED</li> <li>• Shipping weight: 91 lbs.</li> <li>• *This product is available left hinged, product number: ABT-HC-UCFS-0104-LH</li> </ul>