

SPEC SHEET ABT-HC-UCFS-0104-LH

Product Name:	ABT-HC-UCFS-0104-LH
Features:	<ul style="list-style-type: none"> • Premier Countertop Refrigerator Freestanding-LEFT HINGED • 1 Cu. Ft. Capacity • One solid swing door, left hinged • Microprocessor temperature controller • Digital temperature display • Audible and visual high and low temperature alarms • Remote alarm contacts • 2 probes; 1 in air and 1 in glycol bottle • Adjustable operating temperature range: 1°C – 10°C • 2 adjustable shelves • Back wall evaporator cover • Forced draft circulation • Keyed door lock • Integrated door handle • Magnetic door gasket for positive seal • Leveling legs • LED interior light • Probe access port (3/8") • Powder coat exterior • High density urethane foam cabinet and door insulation • HFC-free refrigerant (R600a), foam insulation and packaging • Two year parts and labor warranty, plus an additional three year compressor parts only warranty • Pyxis®, Omnicell® and AcuDose RX® compatible • Exterior dimensions: 17 ¼" W x 19 ¼" D x 21 ¼" H • Note: This unit must have 4" clearance on sides and back for adequate ventilation

SPEC SHEET **ABT-HC-UCFS-0104-LH**

Product Name:	ABT-HC-UCFS-0104-LH
Features:	<ul style="list-style-type: none">• Cycle defrost• 115V, 60 Hz, .90 Amps, 1/8 HP• Energy Consumption 0.71 kWh/day UL/C-UL listed• ENERGY STAR CERTIFIED• Shipping weight: 91 lbs.