

SPEC SHEET **ABT-HC-UCFS-0204G**

Product Name:	ABT-HC-UCFS-0204G
Features:	<ul style="list-style-type: none"> • Premier Undercounter Refrigerator Freestanding • 2.5 Cu. Ft. Capacity • One glass swing door, right hinged, not reversible • Microprocessor temperature controller • Adjustable operating temperature range: 1°C – 10°C • Digital temperature display • Audible and visual high and low temperature alarms • 2 probes; 1 in air and 1 in sample bottle • Remote alarm contacts • Four adjustable shelves (two large, two small) with guard rail on back to prevent product from being pushed against the back wall • Back wall evaporator cover • Forced draft circulation • Keyed door lock • Magnetic door gasket for positive seal • Leveling legs • Probe access port (3/8") • Powder coat exterior • High density urethane foam cabinet and door insulation • HFC-free refrigerant (R600a), foam insulation and packaging • Two year parts and labor warranty • Pyxis®, Omnicell® and AcuDose RX® compatible • Exterior dimensions: 17 ¾" W x 19 ½" D x 28 7/8" H • Note: This unit must have 4" clearance on sides and back for adequate ventilation

SPEC SHEET **ABT-HC-UCFS-0204G**

Product Name:	ABT-HC-UCFS-0204G
Features:	<ul style="list-style-type: none">• Cycle defrost• 115V, 60 Hz, .90 Amps, 1/8 HP• Energy Consumption 0.58 kWh/day• ADA compliant• UL/C-UL listed• ENERGY STAR CERTIFIED• Shipping weight: 97 lbs.