

SPEC SHEET **ABT-HC-UCFS-0420W**

Product Name:	ABT-HC-UCFS-0420W
Features:	<ul style="list-style-type: none"> • Standard Undercounter Freezer Freestanding • 4 Cu. Ft. Capacity • One solid swing door, right hinged, not reversible • Mechanical thermostat • Adjustable operating temperature range: -15°C – -25°C • 2 fixed shelves • Keyed door lock • Magnetic door gasket for positive seal • Recessed door handle • Leveling legs • Probe access port (3/8") • White cabinet and door • High density urethane foam cabinet and door insulation • HFC-free refrigerant (R600a), foam insulation and packaging • One year parts and labor warranty • Exterior dimensions: 20 1/4" W x 21 5/8" D x 31 1/2" H • Note: This unit must have 4" clearance on sides and back for adequate ventilation • Manual defrost • 115V, 60 Hz, .90 Amps, 1/7 HP • UL/C-UL listed • Shipping weight: 94 lbs.