

SPEC SHEET ABT-HC-UCFS-0520

Product Name:	ABT-HC-UCFS-0520
Features:	<ul style="list-style-type: none"> • Premier Undercounter Laboratory Freezer Freestanding • (-20°C) 4 Cu. Ft. Capacity • One solid swing door, right hinged, not reversible • Microprocessor temperature controller; Digital temperature display • Audible and visual high and low temperature alarms; Remote alarm contacts • 2 probes; 1 in air and 1 in sample bottle • Adjustable operating temperature range: -15°C – -25°C • 2 fixed shelves • Bottom footed shelf to allow air circulation • Back wall evaporator cover • Keyed door lock; Magnetic door gasket for positive seal; Recessed door handle • Leveling legs • LED interior light • Probe access port (3/8") • White cabinet and door • High density urethane foam cabinet and door insulation • HFC-free refrigerant (R600a), foam insulation and packaging • Two year parts and labor warranty, plus an additional three year compressor parts only warranty • Exterior dimensions: 23 5/8" W x 24 1/4" D x 32 3/16" H • Note: This unit must have 4" clearance on sides and back for adequate ventilation • Manual defrost • 115V, 60 Hz, 1.4 Amps, 1/8 HP • UL/C-UL listed • Shipping weight: 128 lbs.