

**SPEC SHEET**    **ABT-MFP-30-TS**

<b>Product Name:</b>	ABT-MFP-30-TS
<b>Features:</b>	<ul style="list-style-type: none"> <li>• TempLog Premier Manual Defrost Laboratory Freezer</li> <li>• 30 Cu. Ft. capacity</li> <li>• One solid swing door, right hinged, not reversible</li> <li>• Five inner "fast-freeze compartments" with doors for temperature maintenance and recovery</li> <li>• Microprocessor Temperature Controller</li> <li>• Adjustable Operating Temperature Range: -15°- -25°C</li> <li>• Touch Screen Digital Temperature Display and Alarm Module with battery back-up to support temperature monitoring and alarms</li> <li>• Password-protected Control Panel</li> <li>• Sample simulated digital temperature display; .1° Resolution, adjustable temperature display in Celsius or Fahrenheit</li> <li>• Active Temperature Scale Bar Display</li> <li>• Resettable 1-point Min/Max History</li> <li>• Digital 24-hour Temperature Chart</li> <li>• Temperature Data Logging</li> <li>• USB Port for Data Transfer and Software Updates</li> <li>• Real-time Operational Display Clock</li> <li>• 4-20mA Output</li> <li>• Electronic Troubleshooting Guide</li> <li>• Electronic QR Code to Access Owner's Manual &amp; Technical Support Contact Information</li> <li>• Adjustable Alarm Volume</li> <li>• Audible and Visual High/Low Temperature Alarms</li> <li>• Sensor Error Alarm</li> <li>• Power Failure Alarm</li> </ul>

**SPEC SHEET ABT-MFP-30-TS**

<b>Product Name:</b>	ABT-MFP-30-TS
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Door Ajar Alarm</li> <li>• Low Battery Alarm</li> <li>• Operation Indicator</li> <li>• Alarm Validation (audible, visual, and remote alarm)</li> <li>• Alarm Mute with Ring Back</li> <li>• Remote Alarm Contacts</li> <li>• Four fixed shelves</li> <li>• Keyed door lock</li> <li>• Magnetic door gasket for positive seal</li> <li>• Stainless steel door handle</li> <li>• Caster platform (adds 4" to height)</li> <li>• One probe access port (3/8")</li> <li>• White textured cabinet and door</li> <li>• High density urethane foam cabinet and door insulation</li> <li>• CFC-free refrigerant (R404a), foam insulation and packaging</li> <li>• Three year parts and labor warranty, plus additional two year compressor parts warranty</li> <li>• Exterior dimensions: 35 3/4" W x 34 1/4" D x 76 1/2" H with caster platform; 72 1/2" H without caster platform</li> <li>• Manual defrost</li> <li>• 115V, 60 Hz, 6 Amps, 1/3 HP</li> <li>• UL/C-UL listed</li> <li>• Shipping weight: 462 lbs.</li> </ul>