

SPEC SHEET **ABT-MFP-3030-TS**

Product Name:	ABT-MFP-3030-TS
Features:	<ul style="list-style-type: none"> • TempLog Premier Manual Defrost Laboratory Freezer, -30C • 30 Cu. Ft. capacity • One solid swing door, right hinged, not reversible • Five inner "fast-freeze compartments" with doors for temperature maintenance and recovery • Microprocessor temperature controller: Adjustable operating temperature range: -27°C – -33°C • Touch Screen Digital Temperature Display and Alarm Module with battery back-up to support temperature monitoring and alarms • Password-protected control panel • Sample simulated digital temperature display; .1° Resolution, adjustable temperature display in Celsius or Fahrenheit • Active temperature scale bar display • Resettable 1-point min/max history • Digital 24-hour temperature chart • Temperature data logging • USB port for data transfer and software updates • Real-time operational display clock • 4-20mA output • Electronic troubleshooting guide • Electronic QR Code to access Owner's Manual & technical support contact Information • Adjustable alarm volume • Audible and visual high/low temperature alarms • Sensor error alarm

SPEC SHEET ABT-MFP-3030-TS

Product Name:	ABT-MFP-3030-TS
Features:	<ul style="list-style-type: none"> • Power failure alarm • Door ajar alarm • Low battery alarm • Operation indicator • Alarm validation (audible, visual, and remote alarm) • Alarm mute with ring back • Remote alarm contacts • Four fixed shelves • Keyed door lock • Magnetic door gasket for positive seal • Stainless steel door handle • Caster platform (adds 4" to height) • One probe access port (3/8") • White textured cabinet and door • High density urethane foam cabinet and door insulation • CFC-free refrigerant (R404a), foam insulation and packaging • Three year parts and labor warranty, plus additional two year compressor parts warranty • Exterior dimensions: 35 3/4" W x 34 1/4" D x 76 1/2" H with caster platform; 72 1/2" H without caster platform • Manual defrost • 115V, 60 Hz, 10.7 Amps, 1 HP • UL/C-UL listed • Shipping weight: 462 lbs.