

SPEC SHEET CRT-ABT-HC-S12G

Product Name:	CRT-ABT-HC-S12G
Features:	<ul style="list-style-type: none"> • Upright Controlled Room Temperature Cabinet, Glass Door • 12 Cu. Ft. Capacity • One swing glass door, self-closing (energy efficient, double pane vacuum insulated), right hinged, not reversible • Microprocessor temperature controller • Digital temperature display • Audible and visual high and low temperature alarms • Remote alarm contacts • Adjustable operating temperature range: 68°F–77°F (20°C – 25°C) • Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years certification of calibration, “buffered” probe in the product simulated solution, min/max memory, F/C switchable, field installable, and visual & audible temp alarm • 2 probes; 1 in air and 1 in sample bottle • Three adjustable shelves • Vaccine storage power cord warning label • Forced draft air circulation • Keyed door lock • Magnetic door gasket for positive seal • Heavy duty door handle • Swivel Casters (3”) • LED interior lights are safety shielded and switch controlled • Probe access port (3/8”) • Powder coated steel white interior and exterior • High density urethane foam cabinet insulation • Hydrocarbon, natural refrigerant (R290)

SPEC SHEET CRT-ABT-HC-S12G

Product Name:	CRT-ABT-HC-S12G
Features:	<ul style="list-style-type: none">• One year parts and labor warranty, plus an additional four year compressor parts warranty• Exterior dimensions: 25" W x 29 ¾" D x 65 ¾" H with casters; 62 ¼" H without casters• Cycle defrost• 115V, 60 Hz, 3 Amps, 1/5 HP• ETL listed• Shipping weight: 258 lbs.