

SPEC SHEET PH-ABT-HC-S23G

| | |
|----------------------|---|
| Product Name: | PH-ABT-HC-S23G |
| Features: | <ul style="list-style-type: none"> • Standard Pharmacy/Vaccine Glass Door Refrigerator • 23 Cu. Ft. Capacity • One swing glass door, self-closing (energy efficient, double pane vacuum insulated), right hinged, not reversible • Microprocessor temperature controller • Digital temperature display • Audible and visual high and low temperature alarms • 2 probes; 1 in air and 1 in sample bottle • Remote alarm contacts • Adjustable operating temperature range: 36°F – 46°F (2°C – 8°C) • Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years certification of calibration, “buffered” probe in the product simulated solution, min/max memory. F/C switchable, field installable, and visual & audible temp alarm • Four adjustable shelves (adjustable in 1/2" increments) • Vaccine storage power cord warning label • Forced draft air circulation • Keyed door lock • Magnetic door gasket for positive seal • Heavy duty door handle • Swivel Casters (3 1/2") • LED interior lights are safety shielded and switch controlled • Probe access port (3/8") • Powder coated steel white interior and exterior • High density urethane foam cabinet insulation • Hydrocarbon, natural refrigerant (R290) |

SPEC SHEET PH-ABT-HC-S23G

| | |
|----------------------|--|
| Product Name: | PH-ABT-HC-S23G |
| Features: | <ul style="list-style-type: none">• Pharmacy refrigerator/freezer toolkit and temperature logs• One year parts and labor warranty, plus an additional four year compressor parts warranty• Exterior dimensions: 26 7/8" W x 34 3/4" D x 81 3/4" H with casters; 78 1/4" H without casters• Cycle defrost• 115V, 60 Hz, 3 Amps, 1/5 HP• ETL listed• Shipping weight: 336 lbs. |