

SPEC SHEET PH-ABT-HC-SSP-23

Product Name:	PH-ABT-HC-SSP-23
Features:	<ul style="list-style-type: none"> • Premier Pharmacy/Vaccine Stainless Steel Refrigerator • 23 Cu. Ft. Capacity • One stainless steel swing door, right hinged, not reversible • Microprocessor temperature controller • Digital temperature display • Audible and visual high and low temperature alarms • Remote alarm contacts • 2 probes; 1 in air and 1 in sample bottle • Adjustable operating temperature range: 36°F – 46°F (2°C- 8°C) • Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration, with min/max memory, °F/°C switchable • Three adjustable shelves • Vaccine storage power cord warning label • Forced draft air circulation • Keyed door lock • Magnetic door gasket for positive seal • Integrated door handle • Swivel casters (adds 6 1/2" to height) • Interior light • Probe access port (3/8") • Stainless steel interior floor; Stainless steel exterior • High density urethane foam cabinet insulation • Hydrocarbon, natural refrigerant (R290) • Pharmacy refrigerator/freezer toolkit and temperature logs • Two year parts and labor warranty, plus an additional three year compressor parts warranty

SPEC SHEET PH-ABT-HC-SSP-23

Product Name:	PH-ABT-HC-SSP-23
Features:	<ul style="list-style-type: none">• Exterior dimensions: 26 ¾" W x 34 7/8" D x 83 7/8" H with casters; 77 3/8" H without casters• Cycle defrost• 115V, 60 Hz, 3 Amps, 1/5 HP• ETL listed• Shipping weight: 295 lbs.