

**SPEC SHEET PH-ABT-HC-SSP-72G**

<b>Product Name:</b>	PH-ABT-HC-SSP-72G
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Premier Series Pharmacy Stainless-Steel Glass Door Refrigerator</li> <li>• 72 Cu. Ft. capacity</li> <li>• Three swing glass doors, self-closing, not reversible</li> <li>• 12 ea. adjustable shelves</li> <li>• Digital display microprocessor temperature controller</li> <li>• Adjustable operating temperature: (36°F to 46°F)</li> <li>• Temperature Monitor Device - complies with the current CDC guidelines, 3 years certification of calibration, "buffered" probe in the product simulated solution, min/max memory, °F/°C switchable, field installable, and visual &amp; audible temperature alarms</li> <li>• 2 probes; 1 in air and 1 in sample bottle</li> <li>• Sensor error alarm</li> <li>• Continuous product (simulated) temperature display</li> <li>• Remote alarm contacts</li> <li>• Forced draft air circulation</li> <li>• Keyed door locks</li> <li>• Magnetic door gaskets for positive seal</li> <li>• Integrated door handle</li> <li>• LED interior lights are safety shielded and switch controlled</li> <li>• Probe access port (3/4")</li> <li>• Stainless steel exterior</li> <li>• White-powder coated aluminum interior with 304 stainless steel floor</li> <li>• Easily cleaned smooth interior</li> <li>• Hydrocarbon natural refrigerant (R290)</li> <li>• 2 year parts and labor warranty, plus an additional 3 year compressor parts warranty</li> </ul>

**SPEC SHEET PH-ABT-HC-SSP-72G**

<b>Product Name:</b>	PH-ABT-HC-SSP-72G
<b>Features:</b>	<ul style="list-style-type: none"><li>• 6 swivel casters (5 3/4")</li><li>• Exterior dimensions: 81" W x 35 3/4" D x 83 1/4" H with casters; 77 1/2" H without casters</li><li>• Cycle defrost</li><li>• 115V, 60 Hz, 7.5 Amps, 1/2 HP</li><li>• ETL listed</li><li>• Shipping weight: 786 lbs.</li></ul>