

SPEC SHEET PH-ABT-HC-UCBI-0204-LH

Product Name:	PH-ABT-HC-UCBI-0204-LH
Features:	<ul style="list-style-type: none"> • Premier Pharmacy/Vaccine Undercounter Refrigerator Built-In • 2.5 Cu. Ft. Capacity • One solid swing door, left hinged • Microprocessor temperature controller • Digital temperature display • Adjustable operating temperature range: 36°F – 46°F (2°C – 8°C) • Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration, with min/max memory, °F/°C switchable • Audible and visual high and low temperature alarms • 2 probes; 1 in air and 1 in glycol bottle • Remote alarm contacts • Four adjustable shelves (two large, two small) with guard rail on back to prevent product from being pushed against the back wall • Back wall evaporator cover • Vaccine storage power cord warning label • Forced draft circulation • Self-closing door • Keyed door lock • Magnetic door gasket for positive seal • Leveling legs • LED light • Probe access port (3/8") • Powder coat exterior • High density urethane foam cabinet and door insulation • HFC-free refrigerant (R600a), foam insulation and packaging

SPEC SHEET PH-ABT-HC-UCBI-0204-LH

Product Name:	PH-ABT-HC-UCBI-0204-LH
Features:	<ul style="list-style-type: none"> • Pharmacy refrigerator/freezer toolkit and temperature logs • Two year parts and labor warranty, plus an additional three year compressor parts only warranty • Pyxis®, Omnicell® and AcuDose RX® compatible • Exterior dimensions: 17 3/4" W x 22 3/4" D x 30 3/4" H • Energy Star Certified for Pharmacy Refrigerators and Freezers • Cycle defrost • 115V, 60 Hz, .90 Amps, 1/8 HP • Energy Consumption 0.58 kWh/day • ADA compliant • UL/C-UL listed • Shipping weight: 85 lbs.