

ITEM

75400-714

Capacity	0.52 cu.ft. (14.8 liters)
Interior Dimensions (W x D x H)	12.5" x 6.1" x 6.0" 318 x 154 x 152 mm
Exterior Dimensions (W x D x H)	30.4" x 24.2" x 12.2" 772 x 615 x 310 mm
Dry Weight*	80.8 lb (36.7 kg)
Net Weight**	125.6 lb (57.0 kg)
Shipping Weight	167.0 lb (76.0 kg)
Interior Construction	Aluminum
Exterior Construction	300 Stainless Steel/CR
Wall Clearance, Each Side	4" (102 mm)
Wall Clearance, Rear	4" (102 mm)
Tank Size	14 liters
Tank Fill Capacity	10 liters

ELECTRICAL

Electrical Supply /Hz, AC, 1Ø	50/60 Hz, 100-240 AC, 550 Watts
Plug Supplied	NEMA5-15P, EU1-16P




TEMPERATURE SPECIFICATIONS

Temperature Range	2°C to 70°C
Temperature Uniformity	±0.5°C at 4°C ±0.5°C at 37°C ±1.3°C at 70°C
Temperature Stability - At 37°C	0.1

*Bath empty with no beads loaded.

**Bath filled with beads.

CHG BY	CHG NR	REV	DATE	DESCRIPTION



Lab Armor, LLC www.labarmor.com		PROHIBITARY AND CONFIDENTIAL <small>THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF LAB ARMOR, LLC. IT IS TO BE USED ONLY FOR THE PROJECT AND IN WRITTEN PERMISSION OF LAB ARMOR, LLC. © 2015 LAB ARMOR, LLC.</small>	
DRAWN JAB	NAME JAB	DATE 9/8/2020	MODEL: PELTIER BEAD BATH 14L 100-240V
SCALE: 1:7	SIZE B	DWG. NO. 9021226	REV