

MLATM Precision

VistaLab Catalog #	Color	Volume	Accuracy	Precision
1101	Silver	5µL	±5.00%	3.74%
1021	Orange	10µL	±5.00%	2.60%
1102	Silver	15µL	±5.00%	1.73%
1022	Black	20µL	±3.00%	1.35%
1023	White	25µL	±3.00%	1.12%
1103	Silver	30µL	±3.00%	0.93%
1104	Silver	40µL	±3.00%	0.70%
1024	Green	50µL	±2.00%	0.80%
1106	Silver	60µL	±2.00%	0.67%
1107	Silver	70µL	±2.00%	0.57%
1108	Silver	75µL	±2.00%	0.53%
1109	Silver	80µL	±2.00%	0.67%
1110	Silver	90µL	±2.00%	0.57%
1025	Blue	100µL	±1.00%	0.54%
1111	Silver	120µL	±1.00%	0.36%
1112	Silver	150µL	±1.00%	0.36%
1026	Red	200µL	±1.00%	0.30%
1027	Purple	250µL	±1.00%	0.72%
1028	Brown	300µL	±1.00%	0.60%
1113	Silver	400µL	±1.00%	0.50%
1119	Silver	450µL	±1.00%	0.44%
1029	Yellow	500µL	±1.00%	0.40%
1115	Silver	600µL	±1.00%	1.00%
1505	Silver	700µL	±1.00%	0.86%
1117	Silver	800µL	±1.00%	0.75%
1118	Silver	900µL	±1.00%	0.67%
1030	Gray	1000µL	±1.00%	0.60%

