

# PR1MA™ Orbi-Shaker Large

# PR1MA™



- High capacity, up to 9 x 500 mL, 5 x 1 L or 2 x 2 L
- Horizontally circular 19 mm orbit for aeration & mixing
- Digital speed control up to 300 rpm

Instantly exchange flask clamps with MAGIC Clamp™ accessories.

The powerful PR1MA™ Orbi-Shaker is capable of mixing a 14 lb load capacity up to speeds as high as 300 rpm. The standard platform (included) features a non-slip, rubber-coated surface, ideal for tissue culture flasks, Petri dishes, and staining trays.

The convenient MAGIC Clamp™ platform is also available for use with Erlenmeyer flasks and test tube racks. The unique, magnetic attachment method is the easiest way to instantly change between flask clamps of different sizes.



H1000-MR-CMB



## SPECIFICATIONS

<b>SPEED</b>	Variable, 30 to 300 rpm
<b>MAX. CAPACITY</b>	5 x 1 L   14 lbs   6 kg
<b>DIMENSIONS</b>	11 x 13 x 3 in.
<b>WEIGHT</b>	25 lbs   11.25 kg
<b>ELECTRICAL</b>	120V or 230V, 50-60 Hz 600W
<b>WARRANTY</b>	2-Year Warranty

## ORDERING INFORMATION

ITEM #	DESCRIPTION	QUANTITY
OS-3000	Variable Speed, 13 x 12 in. w/ Rubber Mat Platform	1/ea
OS-PLAT	MAGIC Clamp™ Platform for OS-3000	1/ea
H1000-MR-CMB	MAGIC Clamp™ Combination Clamp	1/ea

**MIDSCI™**  
Laboratory Equipment & Supplies

ORDER TODAY! midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998